

Application Note

Real-time Observation of Neuronal Activity in a Zebrafish Larva



Introduction

Real-time, simultaneous observation of individual neuronal activity over a large area of the brain is required to understand how the brain perceives external sensory information. To achieve this, devices capable of wide-field imaging with high spatial and temporal resolution are necessary.

The combined use of our CSU-W1 confocal scanner and the GCaMP calcium indicator enables imaging of brain activity in zebrafish larvae at a remarkably high resolution.

*Reorganized using the maximum intensity projection (MIP) function.




Fig. 1: Real-time observation of neuronal activity in the optic tectum of zebrafish larva *1



- (a) Zebrafish larva five days after fertilization (□: optic tectum)
- (b) 3D imaging of the brain (tectum) expressing GCaMP
- (c) A GCaMP fluorescence image, single frame extracted from a time-lapse movie.
- (d) Superimposed time-lapse images. Areas with strong calcium reaction are shown in red

Data provided by Professor Koichi Kawakami and Associate Professor Akira Muto, Department of Developmental Genetics, National Institute of Genetics, Japan

Reference: Muto et al., Real-Time Visualization of Neuronal Activity during Perception, Current Biology 23(4):307-311

Experiments

A Zebrafish larva (three to five days after fertilization) expressing the calcium indicator GCaMP in the optic tectum was fixed in agarose, and imaged the neuronal activity in its tectum of both the spontaneous one, and of the response to a visual stimulation.

Sample	UAS (Gal4-UAS system*2): GCaMP7a zebrafish larva three to five days after fertilization	
System		 Zebrafish fixed in agarose
		Confocal unit: CSU-W1 (pinhole: 50μm) Laser: 488nm (solid-state laser) Microscope: Axio Imager (Carl Zeiss) Camera: iXon 888 (Carl Zeiss) Objective lens: W Plan-Apochromat 40x W B-Achroplan 20x Software: Metamorph
Imaging conditions	Continuously imaged the optic tectum (depth: 180μm) for 600-frames at 100msec exposure (total: 152sec., 3.94fps)	

Results

By using the CSU-W1 confocal scanner, a large area of the brain of a zebrafish larva was successfully imaged in a single field, which enabled simultaneous observation of individual neuronal activity within the optic tectum at a high spatial resolution (Fig. 1), and thus, the neurons specifically respond to visual stimuli were identified. Furthermore, it was found that observation of the nerve signals from deep part of the brain (depth of 300μm) are also possible.

Conclusion

Visualization of neural circuit activities in the living brain enables the analysis of functional bonds between neurons. By expressing calcium indicators not only in the optic tectum but also in other brain regions, the CSU-W1 confocal scanner is expected to be a powerful tool for the analysis of functional neural circuits across the different brain regions. Furthermore, the temporal resolution can be increased by selecting a higher speed camera.

*1: Optic Tectum is a major part of the midbrain, and in vertebrates such as fish, etc., it is the main visual processor of the brain.

*2: Gal4-UAS system: A transcription regulation mechanism which can express transgenes in a desired place using the yeast transcription factor Gal4 and its target sequence, UAS(stands for Upstream Activation Sequence).

YOKOGAWA

Bio Solution Center , Life Innovation Bussiness HQ
Yokogawa Electric Corporation

Web site :

<https://www.yokogawa.com/solutions/products-platforms/life-science/>

E-mail: /CSU_livecell_imaging@CS.jp.yokogawa.com

TEL: +81-76-258-7028 FAX: +81-76-258-7029

2-3 Hokuyodai, Kanazawa, Ishikawa 920-0177 Japan

Inquiries