Instruction Manual



Model Y/465-50-SC-S7-□□ pH Electrode for Fermentation Use

1. Specifications

Object of measurement: pH of culture medium in a fermentation

tank

Measurement range : pH0 to pH12

Operating temperature: 0 to 105 °C (when measuring pH)

130°C or less (in sterilization process)

Operating pressure $: 600 \, kPa \, G \, \{6 \, kgf/cm^2 \, G\}$ or less

Temperature sensor : None

Electrode insertion length within a fermentation tank

(leg-length - 50 mm) :

70 mm (leg length: 120 mm) 100 mm(leg length: 150 mm) 150 mm(leg length: 200 mm) 200 mm(leg length: 250 mm)

Converter / transmitter to Be Combined:

PH400G pH Converter, or

PH200G pH Transmitter

Mounting holder : Y/764 Holder for pH electrode for

fermentation use
Y/776 Holder for pH electrode for


fermentation use (removable)
Note: The electrode that can be

incorporated in the Y/776 is limited only to that with an insertion length of 200 mm. (The insertion length within the fermentation tank is

70 mm.)

Applicable electrode cable:

Part number K9148KE : length 1 m Part number K9148KF : length 3 m Part number K9148KG : length 5 m Part number K9148KH : length 10 m




pH Electrode for Fermentation Use and Electrode Cable

2. Model and Codes

Model	Basic Code			Option Code	Description
Y / 465	T				pH electrode for fermentation use
-	- 50				Always - 50
5		·S	C		Always - SC
-		1.	S7		Always - S7
within the			- 07 - 10 - 15 - 20		70 mm 100 mm 150 mm 200 mm

3. Wiring

Connect the cable to the terminals of the PH400G pH converter or the PH200G pH transmitter.



Note: The pH electrode for fermentation use does not incorporate a temperature sensor. Thus, there are no connections to the converter/transmitter terminals T2 and T1. The two cable conductors abbreviated SE are both connected to the SE terminal of the converter/transmitter. (There is no connection to the terminal S.)