

Conductivity standard 15

$\mu\text{S/cm}$

Specification Sheet (Part/REF # 238927)



Stability from 1 to 3 years respectively and conductivity ranges from 1.3 $\mu\text{S/cm}$ to 100 mS/cm

Product Specifications

Accuracy at 25 °C	$\pm 1 \%$
Value at 25 °C	15 $\mu\text{S/cm}$
Package	300 mL glass bottle
Stability Months (after date of manufacturing)	36
Certificate	DFM

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/process-analytics/sensors/238927>

Spec. Version B