


# Paperless Data Acquisition Solution for Pharmaceutical and Bio-Tech Applications

## DAQSTATION DX100/DX200


- ▶ Eliminate paper & pens! DAQSTATION saves critical data to digital storage media.
- ▶ **Batch Option** provides Electronic Signature and Data Security; complies with FDA requirement **21 CFR Part11**
- ▶ Optional time saving **Validation Package**
- ▶ Standard **Ethernet** Interface with **TCP/IP** and **FTP** protocol

### Easy Batch Management




- Easy data management with Batch Number, Lot Number & Comments
- Long duration data storage with 100 MB ZIP
- PC software for data and batch information viewing. Also converts data to Excel for custom reports and printing.


### DX200




### Security with Log in and Password Function



### DX100



### Validation Document for IQ and OQ



### 1. Simple Operation-

- Start/Stop keys control data storage; navigation keys change displays modes.

### 2. Data Security & Accuracy-

- Proprietary binary data file protocol with software security check
- User ID and Password saved with data, prevents unauthorized operation.
- Digital data is far more accurate than a pen trace on paper.


### 3. Low Operating Cost-

- No maintenance- eliminate the cost of charts, pens and ribbon cassettes.
- Eliminate chart storage issues by using DX large capacity storage media.

### 4. Versatile Data Management & Viewing-

- Stand-alone operation; recording & historical viewing on the DX
- Network operation: centralized data management using standard Ethernet interface with FTP TCP/IP, Client, and Server modes
- Easy data viewing and conversion using PC software

# Compliance with FDA Regulations Validation Assistance for Your Process **DAQSTATION DX100/DX200**



## DX can review Historical Data

Standard PC Software for File Viewing, plus conversion to Excel for custom reports and printing.

Optional DAQEXPLORER PC Software for Network File Managing functions using FTP, plus Viewing and Configuration.

DX Batch option code /BT1 for Pharmaceutical & Bio Tech applications  
Please refer to catalogue for details

# YOKOGAWA

**YOKOGAWA ELECTRIC CORPORATION**

Measurement Sales Dept./Phone: 81-422-52-6614, Fax: 81-422-52-6624

Network Solutions Sales Dept./Phone: 81-422-52-6765, Fax: 81-422-52-6793

**YOKOGAWA CORPORATION OF AMERICA**

**YOKOGAWA EUROPE B.V.**

**YOKOGAWA ENGINEERING ASIA PTE. LTD.**

Phone: 1-770-253-7000, Fax: 1-770-251-2088

Phone: 31-33-4622142, Fax: 31-33-4641659

Phone: 65-2419933, Fax: 65-2412606

Subject to change without notice.

[Ed : 02/b] Copyright ©2000

Printed in Japan, 103(YG)

RS-08E