



Model DOX8A Accessory



200ml
POLYETHYLEN BEAKER
Weight : Approx. 100g




ZERO ADJUSTING
REAGENT
Weight : Approx. 600g




POLISHER
Weight : Approx. 100g



ELECTROLYTE
Weight : Approx. 100g



MEMBRANE AS.
Weight : Approx. 100g



SYLINGE
Weight : Approx. 100g

Unless otherwise specified, differences in the dimensions are specified as: General tolerance = \pm (Criteria of tolerance class IT18 in JIS B0401-1998)/2.