


## Dissolved-oxygen Metering System DO30G - NN - 25 - □ □ - FL Dissolved - oxygen Sensor

Unit : mm

≠ Dissolved-oxygen sensor with M4 ring terminals



Model and Suffix Codes	L	Weight (kg)
DO30G - NN - 25 - 03 - FL	3 000	Approx. 0.6
DO30G - NN - 25 - 05 - FL	5 000	Approx. 0.8
DO30G - NN - 25 - 10 - FL	10 000	Approx. 1.4
DO30G - NN - 25 - 15 - FL	15 000	Approx. 2.0
DO30G - NN - 25 - 20 - FL	20 000	Approx. 2.6

Unless otherwise specified, differences in the dimensions are specified as : General tolerance=±(Criteria of tolerance class IT18 in JIS B0401)/2.